



Active Vibration Isolation Element Duo100

- dimensions: 1000x208x125mm
- weight: 32kg
- maximum load capacity: 400kg
- material of housing: aluminum
- surface: powder coated
- configuration: 8 sensor/actuator pairs
- electronics: internal
- setup: minimum 2 elements required, minimum distance 400mm, elements need to be mechanically connected

		Allgemein-toleranzen ISO 2768-f		Werkstück-kanten DIN ISO 13715		Werkstoff:		Maßstab:	
		Datum		Name		Active Vibration Isolation Element Duo100			
Gezeichnet		06.04.2020		U. Stöbener					
Kontrolliert									
Norm						Artikel-Nr.:			
						Zeichnungs-Nr.:			
Status		Änderungen		Datum		Name		F:\CAD\Technische Zeichnungen\Duo\Duo100_Standard.idw	