




D-sub connectors require a space of 100mm



Active Vibration Isolation Platform SW1000-4.0

- dimensions: 1000x800x130mm
- weight: 115kg
- maximum load capacity: 600kg
- material of housing: aluminum
- surface: anodized
- configuration: 4 internal isolation elements (16 sensor/actuator pairs)
- electronics: 4-channel controller
- cables: 4 D-sub cables (3m long)

4-channel controller



		Allgemein-toleranzen ISO 2768-f		Werkstück-kanten DIN ISO 13715		Werkstoff:		Maßstab:	
		Datum		Name		Active Vibration Isolation Platform SW1000-4.0 and Controller			
Gezeichnet		21.04.2010		U. Stöbener					
Kontrolliert									
Norm						Zeichnungs-Nr.:			
Artikel-Nr.:									
Status	Änderungen	Datum	Name	F:\CAD\Technische Zeichnungen\MOD-Sandwich Zeichnungen\Standard 1000x800\SW1000_4_Controller.idw				A3	